



## CARDIOPACS + REPORTING:

Standardized cardiovascular report package conforms to ASE (American Society of Echocardiography) & ICAEL reporting format and guidelines, essential for quality assurance program and accreditation, including:

- General echo report
- Stress echo report
- Cerebrovascular report for carotid
- Vascular Reports
- General Reports

Save time and increase accuracy.

Auto-populate measurements and data from ultrasound studies for timely interpretation.

Requires user log-in and password to operate application.

Single or bidirectional EMR interface via HL7 PDF, text-based HL7, or DSR.

Administrative features include user report audits and controlling application privileges.

Provides secure socket layer connections when accessing the software application from the Internet.

Auto and manual calibration tools include measurement, pan/zoom, annotating, error correction, multiple image / video display and more for identifying regions of interest.

Save reports feature.

Synchronize and shuffle stress echo stages and views.

Report dynamic preview.

LAN, WAN, and Internet enabled for online viewing and report creation.

Interfaces with many existing systems such as Philips, GE, SonoSite, Mindray, Aloka, Biosound, ATL, Siemens, and others.

Save or export images/clip as .JPG, AVI, .BMP files.

Telemedicine capabilities for online, remote viewing of studies and to finalize reports.

Easy to use graphical user-interface with point and click navigation.

Provide referring physicians with studies and reports via:

- VidiStar fax and Autofind function
- Burn CDs and view studies with a standard PC using VidiStar's stand alone viewer

